

# Alcohol-free Surface Disinfectant Wipes

30sec

2in1

Mint

x200

! Always Read Product Label and Material Safety Data Sheet before use. For Professional use only.

- Specially formulated and effective against all enveloped viruses, including Coronavirus, HIV, Hepatitis, Influenza and many other microorganisms. (for full details, see overleaf)
- For surfaces and non-invasive medical devices. Ideal to use on surfaces sensitive to alcohol.
- Wipes have a peppermint aroma and are supplied in packs of 200 wipes.
- Each wipe measures 200 x 200mm.
- 2 in 1 product, wipes contain detergent for added cleaning power.

[g9chemicals.com](http://g9chemicals.com)



Your only solution

# Alcohol-free Surface Disinfectant Wipes

30sec

2in1

Mint

x200

ⓘ Always Read Product Label and Material Safety Data Sheet before use. For Professional use only.

## Contact Times

Bactericidal – 30 sec  
Fungicidal – 30 sec  
Yeasticidal – 30 sec  
Enveloped Viruses – 2 min  
Mycobactericidal – 2 min  
Spores – 5 min

## Disinfectant Test Portfolio

The following test standards have been achieved: EN14885, EN14476, EN14348, EN13727, EN13704, EN13624

### Bactericidal:

Pseudomonas aeruginosa,  
Enterococcus hirae,  
Staphylococcus aureus,  
MRSA (60 sec), VRE (60 sec)

### Mycobactericidal:

Mycobacterium avium

### Yeasticidal:

Candida albicans

### Enveloped Viruses:

Poxviridae, Herpesviridae,  
Filoviridae (e.g. Ebola, Marburg),  
Flavivirus, Hepatitis B virus (HBV),  
Hepatitis C virus (HCV), Hepatitis  
Delta virus (HDV), Influenza virus,  
Paramyxoviridae, Rubella virus,  
Measles virus, Rabies virus,  
Coronavirus (e.g. SARS, MERS),  
Human immunodeficiency virus  
(HIV) & Human T cell leukemia virus  
(HTLV)

### Fungicidal:

Aspergillus niger

### Spores:

Clostridium difficile

## Directions For Use:

Using a new wipe for each application, wipe objects to be disinfected so that they are thoroughly moistened and let the disinfectant take action over a minimum of 60 seconds (to achieve full efficacy, leave for five minutes). Allow surface to dry or wipe with clean tissue.

[g9chemicals.com](http://g9chemicals.com)



G9 Chemicals Ltd, 24 High St,  
Charing, Ashford, Kent,  
TN27 0HX, UK

+44 (0)1580 895 141  
[info@g9chemicals.co.uk](mailto:info@g9chemicals.co.uk)  
[g9chemicals.com](http://g9chemicals.com)