

Instrument Cleaner

1:100

2l

! Always Read Product Label and Material Safety Data Sheet before use. For Professional use only.

- Unique formulation designed to remove proteinaceous and organic material from the surface of instruments and endoscopes.
- Will digest and dissolve blood, mucus, faecal and vomit matter. It safely removes any organic contaminants and proteinaceous soils on all surgical instruments as well as fibre optic endoscopes.
- Dilution Ratio Diluted 1:100 (10ml per litre).
- 2 litre concentrate dilutes to 100 litres of working solution.

Directions For Use:

Manual Cleaning

1. Add 50 ml (1:100) of G9 Instrument Cleaner to 5 litres of water.
2. Submerge instruments to be cleaned as soon as possible after use.
3. Rinse thoroughly in demineralised water before further processing, disinfection or sterilisation

Ultrasonic Cleaning machine

1. Add 50 ml (1:100) of G9 Instrument Cleaner to 5 litres of water in the ultrasonic cleaning tank.
2. Leave instruments in the ultrasonic bath for approximately 15 minutes.

3. Remove instruments from the ultrasonic cleaner and rinse thoroughly with demineralised water before further processing, disinfection or sterilisation.

It is not recommended to exceed soaking times or to leave instruments in the cleaning solution overnight.

Working solution lasts up to 7 days.

g9chemicals.com



Your only solution